

**From:** Board of Pharmacy <pharmacy.subscriberlist@DCA.CA.GOV>  
**Sent:** Tuesday, April 18, 2017 11:00 AM  
**To:** PHARM-GENERAL@DCALISTS.CA.GOV  
**Cc:** Pharmacy\_Subscriberlist@DCA  
**Subject:** Recall notice - B Braun

B Braun is recalling the item/lots below due to the potential for internal leaking at the macro valves of APEX 26-lead transfer set in small percentage (.11%) of product. This recall is to the pharmacy level. Affected product started shipping May 27, 2016.

Please immediately check your inventory, quarantine, discontinue distribution of the affected product, and contact your wholesaler for directions.

Description	Lot #/Exp Date	NDC	UPC
APEX 26-LEAD TRANSFER SET CS4	61492651 01/31/18;	08021002234	46964656101
	61503073 04/30/18;		
	61503074 04/30/18;		
	61503075 04/30/18;		
	61503838 04/30/18;		
	61503079 05/21/18;		
	61503076 05/31/18;		
	61503092 05/31/18;		
	61506054 05/31/18;		
	61503096 06/30/18;		
	61506055 06/30/18;		
	61506056 06/30/18;		
	61506057 06/30/18;		
	61507600 07/31/18;		
	61507601 08/31/18;		
	61507602 08/31/18;		
	61507603		

	08/31/18; 61507612 08/31/18; 61514488 08/31/18; 61507645 12/31/18		
--	---	--	--