
From: Board of Pharmacy <pharmacy.subscriberlist@DCA.CA.GOV>
Sent: Thursday, April 20, 2017 8:23 AM
To: PHARM-GENERAL@DCALISTS.CA.GOV
Cc: Pharmacy_Subscriberlist@DCA
Subject: Recall notice - Abbott Medical Optics Inc.

Abbott Medical Optics Inc. (AMO) is recalling 36 product lots of Healon® OVDs of various models. The Healon® OVD lot number is displayed on the end of each individual unit carton. The OVD lot number is also present on each individual syringe packaging tray and each syringe.

AMO has voluntarily initiated this Action because a remote possibility exists that certain Healon®OVD solutions in the affected lots may contain microscopic glass particles due to damage that occurred at the cylinder neck during the manufacturing process. Use of OVD solutions with glass particles could potentially lead to intraocular injury.

Please immediately check your inventory, quarantine, discontinue distribution of the affected product, and contact your wholesaler for directions.

Model	Lot No.
Healon® 0.55 ml	UB32593
Healon® 0.85 ml	UB32540
Healon® 0.85 ml	UB32602
Healon Duet®	UB32636
Healon Duet®	UB32683
Healon Duet®	UB32719
Healon Duet®	UB32720
Healon Duet®	UB32761
Healon GV® 0.55 ml	UB32596
Healon GV® 0.85 ml	UB32571